

# CERTIFICATE OF ANALYSIS

<b>Lysyl Endopeptidase, MS approved lyophil.</b>	<b>Cat.No. : 20987</b>
	<b>Contr.No.: 180234</b>

Parameter	Method	Specification	Result
Appearance		lyophilisate	corresponds
Activity (AU/vial)		0.03 – 0.07	0.06
Minimum shelf life			22.02.2021
Storage (°C)			-15 to -25

**Unit definition:**

One Amidase Unit (AU) is the amount of enzymes, which will produce 1 µmol of p-Nitroaniline per minute at 30 °C, pH 9.5.

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 08.01.2020**

Christian Monsler

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.