

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

E-64 reinst	Kat.Nr./Cat.No. :	21100
trans-Epoxy succinyl-L-Leucylamido-(4-Guanidino)butane research grade	Lot/Contr.No. :	200493

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₁₅ H ₂₇ N ₅ O ₅	
Molekulargewicht molecular weight		357.4	
Aussehen appearance		weißes krist. Pulver white cryst. powder	entspricht corresponds
Gehalt (%) assay	HPLC	min. 99.0	≥ 99
Mindesthaltbarkeit minimum shelf life			31.05.2022
Lagerung (°C) storage			-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 14.05.2020

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.