

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Fast Red TR-Salz	Kat.Nr./Cat.No. : 21317
Fast Red TR-salt	Lot/Contr.No. : 160516

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₇ H ₆ ClN ₂ ·C ₁₀ H ₇ O ₆ S ₂	
Molekulargewicht molecular weight		440.9	
Aussehen appearance		weißes bis leicht gelbliches Pulver white to light yellow powder	entspricht corresponds
λ_{max} (nm)	0.001 % in H ₂ O	286 – 290	286
A 1cm/ λ_{max}	0.001 % in H ₂ O	min. 0.25	0.562
Wasser (%) water	KF	max. 10.0	6.2
TLC	Kieselgel/silicagel 60 F ₂₅₄ n-Butanol 40 Eisessig/acetic acid 20 Wasser/water 20	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			20.05.2021
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 18.09.2017

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.