

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Fibronectin (human)	Kat.Nr./Cat.No. :	21370
lyophil., analytical grade, sterile	Lot/Contr.No. :	110726

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Molekulargewicht molecular weight		ca. 450 000	
Aussehen appearance		weißes Lyophilisat white lyophilisate	entspricht corresponds
Reinheit (%) purity	SDS-Elektrophorese SDS-electrophoresis	> 90	entspricht corresponds
Endotoxine (ng/mg protein)	LAL		< 4
Sterilität sterility		negativ negative	entspricht corresponds
HBsAg , anti-HIV 1 und 2 HBsAg , anti-HIV 1 and 2		negativ negative	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			10/ 2015
Lagerung (°C) storage			+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 05.02.2015

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.