

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

| | |
|--------------------|--------------------------------|
| Ficoll® 400 | Kat.Nr./Cat.No. : 21373 |
| | Lot/Contr.No. : 150368 |

| Parameter parameter | Methode method | Spezifikation specification | Ergebnis result |
|---|--------------------------|---------------------------------------|---------------------------|
| Molekulargewicht molecular weight | | 400 000 ± 100 000 | 383 856 |
| Aussehen appearance | | weißes Pulver white powder | entspricht corresponds |
| Trocknungsverlust (%) loss on drying | | max. 5.0 | 3.0 |
| Gehalt an NaCl (%) content of NaCl | | max. 1.0 | ≤ 1.0 |
| Mindesthaltbarkeit minimum shelf life | | | 10/ 2019 |
| Lagerung (°C) storage | | | +15 bis/ to +30 |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 01.04.2015

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.