

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

<b>Ficoll® 400</b>	<b>Kat.Nr./Cat.No. :</b>	<b>21373</b>
	<b>Lot/Contr.No. :</b>	<b>180349</b>

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Molekulargewicht</b> molecular weight		400 000 ± 100 000	379 722
<b>Aussehen</b> appearance		weißes Pulver white powder	entspricht corresponds
<b>Trocknungsverlust (%)</b> loss on drying		max. 5.0	3
<b>Gehalt an NaCl (%)</b> content of NaCl		max. 1.0	≤ 1.0
<b>Mindesthaltbarkeit</b> minimum shelf life			28.02.2022
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 28.03.2018**

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.