

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Fluoromount W for microscopy	Kat.Nr./Cat.No. : 21634
Fluoromount W for microscopy	Lot/Contr.No. : 130830

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare, farblose Flüssigkeit clear, colorless liquid	entspricht corresponds
Brechungsindex refractive index	$n_D^{20^\circ C}$		1.3778
Mindesthaltbarkeit minimum shelf life			08/ 2014
Lagerung (°C) storage			+2 bis/ to +8

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 22.07.2013

Dipl.-Ing. (FH) Bernhard Göckel

Jakob Arndt

This report has been computer generated and does not contain a signature.