

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Fluoromount W for microscopy	Kat.Nr./Cat.No. : 21634
Fluoromount W for microscopy	Lot/Contr.No. : 130979

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare, farblose Flüssigkeit clear, colorless liquid	entspricht corresponds
Brechungsindex refractive index	$n_D^{20^\circ C}$		1.3784
Mindesthaltbarkeit minimum shelf life			03/ 2015
Lagerung (°C) storage			+2 bis/ to +8

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 13.08.2014

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

This report has been computer generated and does not contain a signature.