

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Fluoromount W for microscopy	Kat.Nr./Cat.No. : 21634
Fluoromount W for microscopy	Lot/Contr.No. : 140268

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare, farblose Flüssigkeit clear, colorless liquid	entspricht corresponds
Brechungsindex refractive index	$n_D^{20^\circ\text{C}}$		1.3796
Mindesthaltbarkeit minimum shelf life			09/ 2015
Lagerung (°C) storage			+2 bis/ to +8

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 06.02.2015

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.