

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Fluoromount W for microscopy	Kat.Nr./Cat.No. : 21634
Fluoromount W for microscopy	Lot/Contr.No. : 170649

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare, farblose Flüssigkeit clear, colorless liquid	entspricht corresponds
Brechungsindex refractive index	$n_D^{20^\circ\text{C}}$		1.3791
Mindesthaltbarkeit minimum shelf life			31.05.2018
Lagerung (°C) storage			+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 10.07.2017

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.