

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Gelrite®	Kat.Nr./Cat.No. :	22168
	Lot/Contr.No. :	121019

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		weißes Pulver white powder	entspricht corresponds
Gelstärke (g/ cm²) gel strength		min. 700	1180
Trocknungsverlust (%) loss on drying		max. 15	5.7
Transmission (%) Transmittance	0.5 % Gel (20 °C) 497 nm	min. 70	83.9
pH		6.0 – 8.0	7.4
Mindesthaltbarkeit minimum shelf life			06/ 2016
Lagerung (°C) storage			+15 bis/ to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 30.01.2015

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.