

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Glaswatte silyliert reinst	Kat.Nr./Cat.No. :	22367
Glass wool, silanized research grade	Lot/Contr.No. :	120024

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		weiße Fasern white fibres	entspricht corresponds
		hydrophobiert mit Dimethyldichlorsilan hydrophobized with dimethyldichlorosilane	
Anwendung use		zum Säulenabschluß in der Gas- und Flüssig- Chromatographie for plugging columns in gas- and liquid chromatography	
Mindesthaltbarkeit minimum shelf life			01/ 2017
Lagerung (°C) storage			+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 23.01.2012

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.