

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Glaswatte silyliert reinst	Kat.Nr./Cat.No. :	22367
Glass wool, silanized research grade	Lot/Contr.No. :	130758

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		weiße Fasern white fibres	entspricht corresponds
		hydrophobiert mit Dimethyldichlorsilan hydrophobized with dimethyldichlorosilane	
Anwendung use		zum Säulenabschluß in der Gas- und Flüssig- Chromatographie for plugging columns in gas- and liquid chromatography	
Mindesthaltbarkeit minimum shelf life			07/ 2018
Lagerung (°C) storage			+15 bis / to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 03.07.2013

Dipl.-Ing. (FH) Bernhard Göckel

Jakob Arndt

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.