

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Glaswatte silyliert reinst	Kat.Nr./Cat.No. :	22367
Glass wool, silanized research grade	Lot/Contr.No. :	150466

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		weiße Fasern white fibres hydrophobiert mit Dimethyldichlorsilan hydrophobized with dimethyldichlorosilane	entspricht corresponds
Anwendung use		zum Säulenabschluß in der Gas- und Flüssig- Chromatographie for plugging columns in gas- and liquid chromatography	
Mindesthaltbarkeit minimum shelf life			04/ 2020
Lagerung (°C) storage			+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 21.04.2015

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.