

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

<b>Glaswatte silyliert</b> reinst	<b>Kat.Nr./Cat.No. :</b>	<b>22367</b>
Glass wool, silanized research grade	<b>Lot/Contr.No. :</b>	<b>160031</b>

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Aussehen</b> appearance		weiße Fasern white fibres	entspricht corresponds
		hydrophobiert mit Dimethyldichlorsilan hydrophobized with dimethyldichlorosilane	
<b>Anwendung</b> use		zum Säulenabschluß in der Gas- und Flüssig- Chromatographie for plugging columns in gas- and liquid chromatography	
<b>Mindesthaltbarkeit</b> minimum shelf life			12.01.2021
<b>Lagerung (°C)</b> storage			+15 bis / to +30

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b>  <b>Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.</b>          We do not guarantee that the product can be used for a special application.          This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 25.01.2016

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.