

**A N A L Y S E N Z E R T I F I K A T**  
CERTIFICATE OF ANALYSIS

---

<b>Glaswatte, silanisiert, reinst</b>	<b>Kat.Nr./Cat.No. :</b>	<b>22367</b>
<b>Glass wool, silanized, research grade</b>	<b>Lot/Contr.No. :</b>	<b>201017</b>

---

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Aussehen</b> appearance		weiße Fasern white fibres  hydrophobiert durch Silanisierung hydrophobized by silanization	entspricht corresponds
<b>Anwendung</b> use		zum Säulenabschluß in der Gas- und Flüssig-Chromatographie for plugging columns in gas- and liquid-chromatography	
<b>Mindesthaltbarkeit</b> minimum shelf life			29.09.2025
<b>Lagerung (°C)</b> storage			+15 bis / to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

---

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 05.10.2020**

Christian Monsler

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.