

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Glaswatte, silanisiert, reinst	Kat.Nr./Cat.No. :	22367
Glass wool, silanized, research grade	Lot/Contr.No. :	220496

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		weiße Fasern white fibres hydrophobiert durch Silanisierung hydrophobized by silanization	entspricht corresponds
Anwendung use		zum Säulenabschluß in der Gas- und Flüssig-Chromatographie for plugging columns in gas- and liquid-chromatography	
Mindesthaltbarkeit minimum shelf life			12.05.2027
Lagerung (°C) storage			+15 bis / to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 15.06.2022

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.