

CERTIFICATE OF ANALYSIS

Glucose oxidase from <i>Aspergillus niger</i> E.C. 1.1.3.4 lyophil.	Cat.No. : 22739 Contr.No.: 130197
---	--

Parameter	Method	Specification	Result
Molecular weight		ca. 160 000	
Appearance		yellow lyophilisate	corresponds
Solubility	5 mg/ml in 0.1 M potassium phosphate buffer, pH 7.0	clear and yellow	corresponds
Activity (U/mg)	Glucose oxidase	min. 220	294
Minimum shelf life			02/ 2015
Storage (°C)			-15 to -25

Unit definition

1 Unit is that amount of enzymatic activity which catalyzes the oxidation of 1 μ mol glucose to glucuronic acid per minute at 25°C, pH 7.0.

We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.
--

SERVA Electrophoresis GmbH
Quality Control**Printing date:** 18.02.2013

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.