

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Glycerin 87% (w/w) pflanzlich
Glycerol 87% (w/w) from plant

Kat.Nr./Cat.No. : 23177
Lot/Contr.No. : P230155

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₃ H ₈ O ₃	
Molekulargewicht molecular weight		92.09	
Aussehen appearance		farblose, viskose Flüssigkeit colorless, viscous liquid	entspricht corresponds
Brechungsindex refractive index	Ph.Eur.	1.449 – 1.455	1.454
Mindesthaltbarkeit minimum shelf life			14.04.2025
Lagerung (°C) storage			+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 18.04.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.