

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Glyceringelatine nach Kaiser phenolfrei	Kat.Nr./Cat.No. :	23311
Glycerol gelatin after Kaiser phenol-free	Lot/Contr.No. :	130438

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₅ H ₈ O ₂	
Molekulargewicht molecular weight		100.1	
Aussehen appearance		leicht gelblich gefärbte feste Masse slightly yellowish colored firm mass	entspricht corresponds
Brechungsindex refractive index	40°C	1.46 – 1.48	entspricht corresponds
Glyceringehalt (%) glycerol content		ca. 40	entspricht corresponds
Gelatinegehalt (%) gelatin content		ca. 7	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			04/ 2018
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 09.04.2013

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.