

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Glyceringelatine nach Kaiser, phenolfrei	Kat.Nr./Cat.No. :	23311
Glycerol gelatin after Kaiser, phenol-free	Lot/Contr.No. :	140821

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		leicht gelblich gefärbte feste Masse slightly yellowish colored firm mass	entspricht corresponds
Brechungsindex refractive index	40 °C	1.44 – 1.48	entspricht corresponds
Glyceringehalt (%) glycerol content		ca. 40	entspricht corresponds
Gelatinegehalt (%) gelatin content		ca. 7	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			07/ 2019
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 24.07.2014

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.