

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Glyceringelatine nach Kaiser, phenolfrei	Kat.Nr./Cat.No. : 23311
Glycerol gelatin after Kaiser, phenol-free	Lot/Contr.No. : 240198

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		leicht gelblich gefärbte feste Masse slightly yellowish colored firm mass	entspricht corresponds
Brechungsindex refractive index	20 °C	1.41 – 1.47	1.41
Glyceringehalt (%) glycerol content		ca. 40	entspricht corresponds
Gelatinegehalt (%) gelatin content		ca. 7	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			27.02.2029
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 28.02.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.