

CERTIFICATE OF ANALYSIS

Insulin	Cat.No. : 26360
recombinant human	Contr.No.: 130953

Parameter	Method	Specification	Result
Molecular weight		ca. 5800	
Assay (%)	calculated on dry basis	95.0 – 105.0	98.4
Assay (U/mg)	USP calculated on dry basis	min. 27.5	28.5
Assay (U/mg)	as is	N/A	26.3
Loss on drying (%)		max. 10.0	7.8
Bacterial endotoxins (IU/mg)		max. 10	< 0.5
Zinc (%)	calculated on dry basis	0.3 – 0.6	0.4
Minimum shelf life			03/ 2017
Storage (°C)			-15 to -25

The physical/ chemical data corresponds to the requirements of Ph. Eur. and USP.

**We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.**

SERVA Electrophoresis GmbH
Quality Control

Printing date: 20.08.2013

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Steinmann

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.