

CERTIFICATE OF ANALYSIS

Catalase from bovine liver lyophil. salt-free	Cat.No. : 26910 Contr.No.: 120128
---	--

Parameter	Method	Specification	Result
Molecular weight		ca. 240 000	
Appearance		green to light brown lyophilisate	corresponds
Solubility	1 mg/ml 0.05 M potassium- phosphate buffer, pH 7.0	clear and light brown	corresponds
Activity (U/mg)	Catalase	ca. 11 000	11 381
Minimum shelf life			07/ 2013
Storage (°C)			-15 to -25

Unit definition

1 Unit is that amount of enzymatic activity which catalyzes the decomposition of 1 µmol of hydrogen peroxide per minute at 25°C, pH 7.0.

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 09.02.2012

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.