

**CERTIFICATE OF ANALYSIS**

<b>Catalase from bovine liver</b> lyophil. salt-free	<b>Cat.No. : 26910</b> <b>Contr.No.: 121234</b>
---	--

Parameter	Method	Specification	Result
<b>Molecular weight</b>		ca. 240 000	
<b>Appearance</b>		green to light brown lyophilisate	corresponds
<b>Solubility</b>	1 mg/ml 0.05 M potassium- phosphate buffer, pH 7.0	clear and light brown	corresponds
<b>Activity (U/mg)</b>	Catalase	ca. 11 000	11 044
<b>Minimum shelf life</b>			11/ 2014
<b>Storage (°C)</b>			-15 to -25

**Unit definition**

1 Unit is that amount of enzymatic activity which catalyzes the decomposition of 1 µmol of hydrogen peroxide per minute at 25°C, pH 7.0.

<b>We do not guarantee that the product can be used for a special application.</b> <b>This document does not release you from performing the standard control upon receipt of incoming goods.</b>
--

**SERVA Electrophoresis GmbH**  
**Quality Control****Printing date:** 07.11.2012

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.