

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Lecithin aus Eidotter rein	Kat.Nr./Cat.No. :	27608
Lecithin from egg yolk pure	Lot/Contr.No. :	130383

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		gelbbbeige Agglomerate yellowish agglomerate	entspricht corresponds
Phosphatidylcholin (%) phosphatidylcholine		min. 60	> 60
unpolare Lipide (%) unpolar lipids		max. 20.0	< 20
Jodzahl iodine number		min. 60.0	68
Peroxidzahl peroxid value		max. 3	0
Wasser (%) water	KF	max. 2.0	0.3
Mindesthaltbarkeit minimum shelf life			01/ 2015
Lagerung (°C) storage			-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 04.04.2013

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Steinmann

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.