

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Lecithin aus Eidotter rein	Kat.Nr./Cat.No. :	27608
Lecithin from egg yolk pure	Lot/Contr.No. :	170130

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		gelbb beige Agglomerate yellowish agglomerate	entspricht corresponds
Phosphatidylcholin (%) phosphatidylcholine		min. 70	70
Jodzahl iodine number		63 – 69	68
Peroxidzahl peroxid value		max. 5	0
Wasser (%) water	KF	max. 2.0	0.4
Mindesthaltbarkeit minimum shelf life			31.07.2019
Lagerung (°C) storage			-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 06.02.2017

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.