

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Lecithin aus Eidotter rein Lecithin from egg yolk pure	Kat.Nr./Cat.No. : 27608 Lot/Contr.No. : 180472
--	---

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		gelbbeige Agglomerate yellowish agglomerate	entspricht corresponds
Phosphatidylcholin (%) phosphatidylcholine		min. 70	70
Jodzahl iodine number		63 – 69	67
Peroxidzahl peroxid value		max. 5	0
Wasser (%) water	KF	max. 2.0	0.5
Mindesthaltbarkeit minimum shelf life			31.10.2020
Lagerung (°C) storage			-15 bis/ to -25

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 30.04.2018

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.