

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

Lectin aus <i>Canavalia ensiformis</i> lyophil.	Kat.Nr./Cat.No. : 27648
Lectin from <i>Canavalia ensiformis</i> lyophil.	Lot/Contr.No. : 220333

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Molekulargewicht molecular weight		110 000	
Aussehen appearance		weißes Lyophilisat white lyophilisate	entspricht corresponds
Löslichkeit solubility	1 mg/ ml PBS	leicht trübe, farblose Lösung slightly cloudy, colorless solution	entspricht corresponds
Protein (%)	aus/ from A _{280 nm}	min. 93	98
Elektrophorese electrophoresis	SDS-PAGE	4 Untereinheiten 4 subunits	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			30.04.2026
Lagerung (°C) storage			+2 bis/ to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 14.05.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.