

ANALYSENZERTIFIKAT
 CERTIFICATE OF ANALYSIS

3-Methacryloxypropyltrimethoxysilan	Kat.Nr./Cat.No. :	28739
3-Methacryloxypropyltrimethoxysilane (Bind-Silane)	Lot/Contr.No. :	190471

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₁₀ H ₂₀ O ₅ Si	
Molekulargewicht molecular weight		248.4	
Aussehen appearance		klare, farblose Flüssigkeit clear, colorless liquid	entspricht corresponds
Gehalt (%) assay	GC	98.0 – 100.0	99.1
Brechungsindex refractive index	20 °C	1.430 – 1.433	1.4315
Farbzahl color	APHA	max. 30	14
Viskosität (mm²/s) viscosity		2.3 – 2.5	2.4
Mindesthaltbarkeit minimum shelf life			07.05.2021
Lagerung (°C) storage			+2 bis/ to +8

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control
Druckdatum/ printing date: 09.05.2019

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.