

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

N,N'-Methylenbisacrylamid 2X p.A.	Kat.Nr./Cat.No. : 29195
Bis(acrylamido)methane 2X analytical grade	Lot/Contr.No. : 221201

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₇ H ₁₀ N ₂ O ₂	
Molekulargewicht molecular weight		154.2	
Aussehen appearance		weiße Nadeln white needles	entspricht corresponds
Gehalt (%) assay	HPLC	min. 98.0	99.1
A 1 cm/ 290 nm	1% in H ₂ O	max. 0.25	0.172
pH	1% in H ₂ O	5.5 – 7.5	5.6
Löslichkeit solubility	5% in H ₂ O/ 50 °C	klar, farblos clear, colorless	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			07.12.2026
Lagerung (°C) storage			+2 bis/ to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 13.12.2022

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.