

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Natriumazid reinst	Kat.Nr./Cat.No. :	30175
Sodium azide research grade	Lot/Contr.No. :	150564

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		NaN ₃	
Molekulargewicht molecular weight		65.01	
Aussehen appearance		weiße Kristalle white crystals	entspricht corresponds
Gehalt (%) assay	Titration	min. 99.0	99.5
Trocknungsverlust loss on drying	105 °C, 4 h	max. 0.1	0.07
Mindesthaltbarkeit minimum shelf life			05/ 2020
Lagerung (°C) storage			+15 bis/ to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 08.05.2015

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.