

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Natriumazid reinst	Kat.Nr./Cat.No. :	30175
Sodium azide research grade	Lot/Contr.No. :	200822

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		NaN ₃	
Molekulargewicht molecular weight		65.01	
Aussehen appearance		weiße Kristalle white crystals	entspricht corresponds
Gehalt (%) assay	Titration	min. 99.0	100.1
Trocknungsverlust loss on drying	105 °C, 4 h	max. 0.1	0.06
Mindesthaltbarkeit minimum shelf life			10.08.2025
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 11.08.2020

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.