

ANALYSENZERTIFIKAT

CERTIFICATE OF ANALYSIS

Nonenylbernsteinsäurenanhydrid rein	Kat.Nr./Cat.No. :	30812
Nonenylsuccinic anhydride pure	Lot/Contr.No. :	071018

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₁₃ H ₂₀ O ₃	
Molekulargewicht molecular weight		224.3	
Aussehen appearance		klare, schwach gelbe, viskose Flüssigkeit clear, slightly yellow, viscous fluid	entspricht corresponds
Farbe color value	Gardner Skala gardner scale	max. 5	< 5
Säurezahl acid number		470 – 495	494
Säuregehalt (%) acidity	als/ as nonenyl succinic acid	max. 2	< 2
Identität identity	IR	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			26.11.2016
Lagerung (°C) storage			+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 28.08.2015

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.