

# ANALYSENZERTIFIKAT

## CERTIFICATE OF ANALYSIS

|                 |                   |        |
|-----------------|-------------------|--------|
| Norit® A pract. | Kat.Nr./Cat.No. : | 30890  |
| Norit® A pract. | Lot/Contr.No. :   | 130012 |

| Parameter<br>parameter                                                            | Spezifikation<br>specification                                                         | Ergebnis<br>result        |
|-----------------------------------------------------------------------------------|----------------------------------------------------------------------------------------|---------------------------|
| <b>Beschreibung</b><br>description                                                | Aktivkohle aus Torf,<br>säuregewaschen<br>activated charcoal from peat,<br>acid washed | entspricht<br>corresponds |
| <b>Aussehen</b><br>appearance                                                     | feines schwarzes Pulver<br>fine black powder                                           | entspricht<br>corresponds |
| <b>Wasser (%)</b><br>water                                                        | max. 10.0                                                                              | 4.0                       |
| <b>Asche (%)</b><br>ash                                                           |                                                                                        | 5.0                       |
| <b>Salzsäurelösliche Bestandteile (%)</b><br>hydrochloric acid soluble components | max. 1.0                                                                               | 0.2                       |
| <b>Methylenblau-Adsorption (g/100g)</b><br>methylene blue adsorption              |                                                                                        | 12                        |
| <b>Molassezahl</b><br>molasses number                                             | max. 380                                                                               | 320                       |
| <b>Mindesthaltbarkeit</b><br>minimum shelf life                                   |                                                                                        | 01/ 2018                  |
| <b>Lagerung (°C)</b><br>storage                                                   |                                                                                        | +15 bis/ to +30           |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.  
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 07.01.2013

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**

This report has been computer generated and does not contain a signature.