

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Nycodenz® 60% (w/v) Lösung in Wasser
solution in water

Kat.Nr./Cat.No. : 31000
Lot/Contr.No. : P250227

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		$C_{19}H_{26}I_3N_3O_9$	
Molekulargewicht molecular weight		821.1	
Aussehen appearance		klare, farblose Lösung clear, colorless solution	entspricht corresponds
Dichte (g/ml) density	20 °C	1.310 ± 0.002	entspricht corresponds
Brechungsindex refractive index	20 °C	1.4273 ± 0.0003	1.4270
Mindesthaltbarkeit minimum shelf life			16.07.2027
Lagerung (°C) storage			+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 17.07.2025

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.