

**CERTIFICATE OF ANALYSIS**

---

<b>Pepsin porcine</b>	<b>Cat.No. : 31820</b>
E.C. 3.4.23.1. 2x cryst. lyophil.	<b>Contr.No.: 130022</b>

---

<b>Parameter</b>	<b>Method</b>	<b>Specification</b>	<b>Result</b>
<b>Molecular weight</b>		ca. 36 000	
<b>Appearance</b>		white lyophilisate	corresponds
<b>Solubility</b>	1 mg/ml 0.01 N HCl	clear and colorless	corresponds
<b>Activity (milliAnson-U/mg)</b>	Pepsin by Anson	ca. 15	12.5
<b>Activity (PU/mg)</b>	„Pepsin Unit“ according to Anson	ca. 0.015	0.0125
<b>Minimum shelf life</b>			01/ 2014
<b>Storage (°C)</b>			-15 to -25

**Unit definition**

1 milliAnson-Unit is equivalent to 1µmol of Folin-positive amino acids calculated as tyrosine, released from denaturated hemoglobin per minute at 37 °C, pH 2.0.

<b>We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</b>
--

---

**SERVA Electrophoresis GmbH**  
**Quality Control****Printing date: 15.01.2013**

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.