

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Phenylmethansulfonylfluorid reinst	Kat.Nr./Cat.No. :	32395
Phenylmethylsulfonyl fluoride research grade	Lot/Contr.No. :	130539

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₇ H ₇ FO ₂ S	
Molekulargewicht molecular weight		174.19	
Aussehen appearance		weiße Kristalle white crystals	entspricht corresponds
Gehalt (%) assay	GC	min. 99.0	99.5
Schmelzpunkt (°C) melting point		90 – 94	92 – 93
Identität identity	IR	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			05/ 2018
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 03.05.2013

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Steinmann

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.