

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

Phenylmethansulfonylfluorid	Kat.Nr./Cat.No. : 32395.54
Phenylmethylsulfonyl fluoride	Lot/Contr.No. : 180921

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₇ H ₇ FO ₂ S	
Molekulargewicht molecular weight		174.19	
Aussehen appearance		weiße Kristalle white crystals	entspricht corresponds
Gehalt (%) assay	GC	min. 99.0	99.2
Schmelzpunkt (°C) melting point		90 – 94	93
Identität identity	IR	entspricht corresponds	entspricht corresponds
Herstelldatum manufacturing date			27.08.2018
Mindesthaltbarkeit minimum shelf life			27.08.2023
Lagerung (°C) storage			+15 bis/ to +30

Country of origin: Germany

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 30.08.2018

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.