

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Phenylmethansulfonylfluorid	Kat.Nr./Cat.No. :	32395
Phenylmethylsulfonyl fluoride	Lot/Contr.No. :	230616

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₇ H ₇ FO ₂ S	
Molekulargewicht molecular weight		174.19	
Aussehen appearance		weiße Kristalle white crystals	entspricht corresponds
Gehalt (%) assay	GC	min. 99.0	99.7
Schmelzpunkt (°C) melting point		90 – 94	92 – 93
Identität identity	IR	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			22.05.2028
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 20.07.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.