

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Polyamid 6 Pulver reinst	Kat.Nr./Cat.No. :	33143
Nylon 6, research grade	Lot/Contr.No. :	131338

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		weißes Pulver white powder	entspricht corresponds
Korngröße (mm) particle size		0.05 – 0.16	entspricht corresponds
pH	10 % in H ₂ O	6.0 – 7.0	6.7
Mindesthaltbarkeit minimum shelf life			12/ 2018
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 02.12.2013

Dipl.-Ing. (FH) Bernhard Göckel

Jakob Arndt

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.