

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Polyvinylpyrrolidon 15 pract.	Kat.Nr./Cat.No. :	33422
Polyvinylpyrrolidone 15	Lot/Contr.No. :	130964

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Molekulargewicht molecular weight		ca. 10 000	
Aussehen appearance		weißes Pulver white powder	entspricht corresponds
Farbe color	Hazen 5 % Lösung/ solution	max. 100	30 – 40
Wasser (%) water	KF	max. 5.0	2.4
pH	5 % in H ₂ O	3.0 – 5.0	4.1
Asche (%) residue on ignition		max. 0.02	0.00
Mindesthaltbarkeit minimum shelf life			26.04.2017
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 13.04.2016

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.