

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

Polyvinylpyrrolidon 15 pract.	Kat.Nr./Cat.No. : 33422
Polyvinylpyrrolidone 15	Lot/Contr.No. : 211287

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Molekulargewicht molecular weight		ca. 10 000	
Aussehen appearance		weißes Pulver white powder	entspricht corresponds
Farbe color	Hazen 5% Lösung/ solution	max. 100	30
Wasser (%) water	KF	max. 5.0	2.0
pH	5% in H ₂ O	3.0 – 5.0	3.9
Asche (%) residue on ignition		max. 0.02	0.01
Mindesthaltbarkeit minimum shelf life			31.07.2025
Lagerung (°C) storage			+15 bis/ to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control
Druckdatum/ printing date: 29.06.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.