

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Polyvinylsulfat·K-salz rein	Kat.Nr./Cat.No. : 33426
Polyvinylsulphate·K-salt pure	Lot/Contr.No. : 150489

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		$(C_2H_3O_4SK)_n$ $n \geq 1500$	
Molekulargewicht molecular weight		$(162.21)_n$	
Aussehen appearance		elfenbeinfarbenes bis bräunliches Pulver ivory-colored to brownish powder	entspricht corresponds
Veresterungsgrad (%) esterification degree	nach Trocknung after drying	min. 90.0	94.6
Löslichkeit solubility	1.2 mg/ ml in H ₂ O	klare, farblose bis leicht gelbliche Lösung clear, colourless to pale yellow solution	entspricht corresponds
Identität identity	IR	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			04/ 2017
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 17.04.2015

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.