

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Polyvinylsulfat-K-salz rein	Kat.Nr./Cat.No. :	33426
Polyvinylsulphate-K-salt pure	Lot/Contr.No. :	160701

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		$(C_2H_3O_4SK)_n$ $n \geq 1\ 500$	
Molekulargewicht molecular weight		$(162.21)_n$	
Aussehen appearance		elfenbeinfarbenes bis bräunliches Pulver ivory-colored to brownish powder	entspricht corresponds
Veresterungsgrad (%) esterification degree	nach Trocknung after drying	min. 90.0	96.0
Löslichkeit solubility	1.2 mg/ ml in H ₂ O	klare, farblose bis leicht gelbliche Lösung clear, colourless to pale yellow solution	entspricht corresponds
Identität identity	IR	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			14.07.2018
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 19.07.2016

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.