

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

Polyvinylsulfat-K-salz rein	Kat.Nr./Cat.No. :	33426
Polyvinylsulphate-K-salt pure	Lot/Contr.No. :	190591

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		weißes bis braunes Pulver white to brown powder	entspricht corresponds
Löslichkeit* solubility	1.2 mg/ ml in H ₂ O	klare bis schwach trübe, farblose bis leicht gelbliche Lösung clear to slightly hazy, colourless to pale yellow solution	entspricht corresponds
Kalium (%) potassium	K	min. 21.0	25.3
Identität identity	IR	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			12.06.2022
Lagerung (°C) storage			+15 bis/ to +30

* **Die Lösung kann vereinzelt unlösliche Partikel enthalten.**
 The solution may occasionally contain insoluble particles.

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 07.05.2020

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.