

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

| | |
|---|--------------------------------|
| Ponceau S Lösung für Elektrophorese (0.2%) | Kat.Nr./Cat.No. : 33427 |
| Ponseau S solution for electrophoresis (0.2%) | Lot/Contr.No. : P150459 |

| Parameter parameter | Methode method | Spezifikation specification | Ergebnis result |
|---|--|--------------------------------|---------------------------|
| Aussehen appearance | | rote Lösung red solution | entspricht corresponds |
| Konzentration concentration | | 0.2 % in 3 % TCA | entspricht corresponds |
| A 1 cm/ 524 nm | 1:100 verd. in H ₂ O 1:100 diluted in H ₂ O | min. 0.9 | 1.048 |
| Proteinfärbung protein staining | CAF | entspricht corresponds | entspricht corresponds |
| Mindesthaltbarkeit minimum shelf life | | | 04.12.2018 |
| Lagerung (°C) storage | | | +15 bis/ to +30 |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 04.12.2015

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.