

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

<b>Ponceau S Lösung für Elektrophorese (0.2%)</b>	<b>Kat.Nr./Cat.No. : 33427</b>
Ponseau S solution for electrophoresis (0.2%)	<b>Lot/Contr.No. : P170073</b>

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance		rote Lösung red solution	entspricht corresponds
<b>Konzentration</b> concentration		0.2% in 3% TCA	entspricht corresponds
<b>A 1 cm/ 524 nm</b>	1:100 verd. in H <sub>2</sub> O 1:100 diluted in H <sub>2</sub> O	min. 0.9	1.058
<b>Proteinfärbung</b> protein staining	CAF	entspricht corresponds	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life			13.02.2020
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 20.02.2017**

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.