

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

Ponceau S Lösung für Elektrophorese (0.2%)	Kat.Nr./Cat.No. : 33427
Ponseau S solution for electrophoresis (0.2%)	Lot/Contr.No. : P180160

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		rote Lösung red solution	entspricht corresponds
Konzentration concentration		0.2% in 3% TCA	entspricht corresponds
A 1 cm/ 524 nm	1:100 verd. in H ₂ O 1:100 diluted in H ₂ O	min. 0.9	0.998
Proteinfärbung protein staining	CAF	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			24.04.2021
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 03.05.2018

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.