

**A N A L Y S E N Z E R T I F I K A T**  
 CERTIFICATE OF ANALYSIS

<b>Ponceau S Lösung für Elektrophorese (0.2%)</b>	<b>Kat.Nr./Cat.No. : 33427</b>
Ponseau S solution for electrophoresis (0.2%)	<b>Lot/Contr.No. : P180235</b>

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
<b>Aussehen</b> appearance		rote Lösung red solution	entspricht corresponds
<b>Konzentration</b> concentration		0.2% in 3% TCA	entspricht corresponds
<b>A 1 cm/ 524 nm</b>	1:100 verd. in H <sub>2</sub> O 1:100 diluted in H <sub>2</sub> O	min. 0.9	1.006
<b>Proteinfärbung</b> protein staining	CAF	entspricht corresponds	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life			03.07.2021
<b>Lagerung (°C)</b> storage			+15 bis/ to +30

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b>  <b>Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.</b>          We do not guarantee that the product can be used for a special application.          This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 04.07.2018**

Christian Monsler

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.